











Product		Cat#	Units/ Case*	Blood Pressure	Skin Pressure	Circumference
SMALL Pink		PRH-028-PI-01A	10	< 130 mmHg	197±24 mmHg	14-28 cm
MEDIUM Yellow		PRH-040-YE-01A	10	< 190 mmHg	254±43 mmHg	24-40 cm
MEDIUM Red		PRH-040-RE-01A	10	< 160 mmHg	246±42 mmHg	24-40 cm
MEDIUM Green		PRH-040-GR-01A	10	< 130 mmHg	234±42 mmHg	24-40 cm
LARGE Brown		PRH-060-BR-01A	10	< 190 mmHg	364±49 mmHg	30-55 cm
LARGE Orange		PRH-060-OR-01A	10	< 160 mmHg	329±47 mmHg	30-55 cm
LARGE Blue		PRH-060-BL-01A	10	< 130 mmHg	309±40 mmHg	30-55 cm
X-LARGE- Black&White		PRH-090-BW-01A	10	< 160 mmHg	320±19 mmHg	50-85 cm
MODEL F		PRH-035-FA-01A	25	< 160 mmHg	248±33 mmHg	14-34 cm
MODEL A		PRH-032-MA-01A	12	< 160 mmHg	176±9 mmHg	22-32 cm

*HemaClear® is sold per case